

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\readout_time\no_PF_no_iPAT
 TA:0:41 PAT:Off Voxel size:3.0×3.0×4.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	25
Dist. factor	25 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2800 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HE1-4

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	12
Delay in TR	1000 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	25
Dist. factor	25 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HE1	On
HE3	On
NE1	Off
HE2	On
HE4	On
NE2	Off
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	124 mm
Frequency 1H	123.256091 MHz
Correction factor	1
SincRFPulse 1H	273.137 V

	Gain	High
	Table position	0 mm
	Img. Scale. Cor.	1.000
Physio	1st Signal/Mode	None
Inline	Inline Composing	Off
	Distortion correction	Off
Sequence	Introduction	On
	Averaging mode	Long term
	Multi-slice mode	Interleaved
	Bandwidth	2442 Hz/Px
	Free echo spacing	Off
	Echo spacing	0.49 ms
	EPI factor	64
	RF pulse type	Normal
	Gradient mode	Fast
	Excitation	Standard
	TX/RX delta frequency	0 Hz
	TX Nucleus	None
	TX delta frequency	0 Hz
	Coil elements	HE1-4
	Acquisition duration	0 ms
BOLD	Delay in TR	1000 ms
	GBP	Off
	PBP	Off
	TTP	Off
	Original images	On
	relCBV	Off
	relCBF	Off
	relMTT	Off
	Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\readout_time\6_8_PF_no_iPAT
 TA:0:41 PAT:Off Voxel size:3.0×3.0×4.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	25
Dist. factor	25 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2800 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HE1-4

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	12
Delay in TR	1000 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	25
Dist. factor	25 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HE1	On
HE3	On
NE1	Off
HE2	On
HE4	On
NE2	Off
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	124 mm
Frequency 1H	123.256091 MHz
Correction factor	1
SincRFPulse 1H	273.137 V

	Gain	High
	Table position	0 mm
	Img. Scale. Cor.	1.000
Physio	1st Signal/Mode	None
Inline	Inline Composing	Off
	Distortion correction	Off
Sequence	Introduction	On
	Averaging mode	Long term
	Multi-slice mode	Interleaved
	Bandwidth	2442 Hz/Px
	Free echo spacing	Off
	Echo spacing	0.49 ms
	EPI factor	64
	RF pulse type	Normal
	Gradient mode	Fast
	Excitation	Standard
	TX/RX delta frequency	0 Hz
	TX Nucleus	None
	TX delta frequency	0 Hz
	Coil elements	HE1-4
	Acquisition duration	0 ms
BOLD	Delay in TR	1000 ms
	GBP	Off
	PBP	Off
	TTP	Off
	Original images	On
	relCBV	Off
	relCBF	Off
	relMTT	Off
	Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\\Research\\Rorden_Prisma\\readout_time\\no_PF_2_iPAT
TA:0:44 PAT:2 Voxel size:3.0×3.0×4.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	Off
Start measurements	single

Routine

Nr. of slice groups	1
Slices	25
Dist. factor	25 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2800 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HE1-4

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	12
Delay in TR	1000 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Separate
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	25
Dist. factor	25 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HE1	On
HE3	On
NE1	Off
HE2	On
HE4	On
NE2	Off
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	124 mm
Frequency 1H	123.256091 MHz
Correction factor	1
SincRFPulse 1H	273.137 V

	Gain	High
	Table position	0 mm
	Img. Scale. Cor.	1.000
Physio	1st Signal/Mode	None
Inline	Inline Composing	Off
	Distortion correction	Off
Sequence	Introduction	On
	Averaging mode	Long term
	Multi-slice mode	Interleaved
	Bandwidth	2442 Hz/Px
	Free echo spacing	Off
	Echo spacing	0.49 ms
	EPI factor	64
	RF pulse type	Normal
	Gradient mode	Fast
	Excitation	Standard
	TX/RX delta frequency	0 Hz
	TX Nucleus	None
	TX delta frequency	0 Hz
	Coil elements	HE1-4
	Acquisition duration	0 ms
BOLD	Delay in TR	1000 ms
	GBP	Off
	PBP	Off
	TTP	Off
	Original images	On
	relCBV	Off
	relCBF	Off
	relMTT	Off
	Distortion Corr.	Off

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\readout_time\no_PF_no_iPAT
 TA:0:16 PAT:Off Voxel size:3.4×3.4×2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	50
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	216 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	3724 ms
TE	35.00 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HE1-4

Contrast

MTC	Off
Magn. preparation	None
Flip angle	72 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	50
Dist. factor	20 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HE1	On
HE3	On
NE1	Off
HE2	On
HE4	On
NE2	Off
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	On - Coil Memory
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	216 mm
A >> P	216 mm
F >> H	120 mm
Frequency 1H	123.256091 MHz
Correction factor	1
SincExcRF 1H	84.755 V

	Gain	High
	Table position	0 mm
	Img. Scale. Cor.	1.000
Physio	1st Signal/Mode	None
	Magn. preparation	None
Inline	Inline Composing	Off
	Distortion correction	Off
Sequence	Introduction	On
	Averaging mode	Long term
	Multi-slice mode	Interleaved
	Contrasts	1
	Bandwidth	2170 Hz/Px
	Flow comp.	No
	Free echo spacing	Off
	Echo spacing	0.53 ms
	EPI factor	64
	Gradient mode	Performance
	RF spoiling	Off
	Physio recording	Off
	Triggering scheme	Standard
	TX/RX delta frequency	0 Hz
	TX Nucleus	None
	TX delta frequency	0 Hz
	Coil elements	HE1-4
	Acquisition duration	0 ms
BOLD	GLM Statistics	Off
	Dynamic t-maps	Off
	Ignore meas. at start	0
	Ignore after transition	0
	Model transition states	Off
	Temp. highpass filter	Off
	Threshold	4.00
	Paradigm size	3
	Motion correction	Off
	Spatial filter	Off
	Delay in TR	0 ms
	Distortion Corr.	Off
	Contrasts	1

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\readout_time\68_PF_no_iPAT
 TA:0:16 PAT:Off Voxel size:3.4×3.4×2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	50
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	216 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	3724 ms
TE	35.00 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HE1-4

Contrast

MTC	Off
Magn. preparation	None
Flip angle	72 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	50
Dist. factor	20 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HE1	On
HE3	On
NE1	Off
HE2	On
HE4	On
NE2	Off
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	On - Coil Memory
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	216 mm
A >> P	216 mm
F >> H	120 mm
Frequency 1H	123.256091 MHz
Correction factor	1
SincExcRF 1H	84.755 V

	Gain	High
	Table position	0 mm
	Img. Scale. Cor.	1.000
Physio	1st Signal/Mode	None
	Magn. preparation	None
Inline	Inline Composing	Off
	Distortion correction	Off
Sequence	Introduction	On
	Averaging mode	Long term
	Multi-slice mode	Interleaved
	Contrasts	1
	Bandwidth	2170 Hz/Px
	Flow comp.	No
	Free echo spacing	Off
	Echo spacing	0.53 ms
	EPI factor	64
	Gradient mode	Performance
	RF spoiling	Off
	Physio recording	Off
	Triggering scheme	Standard
	TX/RX delta frequency	0 Hz
	TX Nucleus	None
	TX delta frequency	0 Hz
	Coil elements	HE1-4
	Acquisition duration	0 ms
BOLD	GLM Statistics	Off
	Dynamic t-maps	Off
	Ignore meas. at start	0
	Ignore after transition	0
	Model transition states	Off
	Temp. highpass filter	Off
	Threshold	4.00
	Paradigm size	3
	Motion correction	Off
	Spatial filter	Off
	Delay in TR	0 ms
	Distortion Corr.	Off
	Contrasts	1

SIEMENS MAGNETOM Prisma_fit syngo MR D13D

\\USER\Research\Rorden_Prisma\readout_time\no_PF_2_iPAT
 TA:0:24 PAT:2 Voxel size:3.4×3.4×2.0 mm Rel. SNR:1.00 :epfid

Properties

Prio Recon	Off
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Wait for user to start	On
Start measurements	single

Routine

Nr. of slice groups	1
Slices	50
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	216 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	3724 ms
TE	35.00 ms
Multi-band accel. factor	1
Filter	None
Coil elements	HE1-4

Contrast

MTC	Off
Magn. preparation	None
Flip angle	72 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Measurements	3
Delay in TR	0 ms
Reconstruction	Magnitude
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Single-shot
Distortion Corr.	Off
Hamming	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Nr. of slice groups	1
Slices	50
Dist. factor	20 %
Position	Isocenter
Phase enc. dir.	A >> P
Phase oversampling	0 %
Multi-slice mode	Interleaved
Series	Interleaved
Nr. of sat. regions	0
Position mode	L-P-H
Fat suppr.	Fat sat.
Special sat.	None
Special sat.	None
Set-n-Go Protocol	Off
Table position	P
Inline Composing	Off

System

Body	Off
HE1	On
HE3	On
NE1	Off
HE2	On
HE4	On
NE2	Off
SP5	Off
SP6	Off
SP7	Off
SP8	Off
SP1	Off
SP2	Off
SP3	Off
SP4	Off
Position mode	L-P-H
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Coil Select Mode	On - Coil Memory
B0 Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto
? Ref. amplitude 1H	0.000 V
Position	Isocenter
Rotation	0.00 deg
R >> L	216 mm
A >> P	216 mm
F >> H	120 mm
Frequency 1H	123.256091 MHz
Correction factor	1
SincExcRF 1H	84.755 V

	Gain	High
	Table position	0 mm
	Img. Scale. Cor.	1.000
Physio	1st Signal/Mode	None
	Magn. preparation	None
Inline	Inline Composing	Off
	Distortion correction	Off
Sequence	Introduction	On
	Averaging mode	Long term
	Multi-slice mode	Interleaved
	Contrasts	1
	Bandwidth	2170 Hz/Px
	Flow comp.	No
	Free echo spacing	Off
	Echo spacing	0.55 ms
	EPI factor	64
	Gradient mode	Performance
	RF spoiling	Off
	Physio recording	Off
	Triggering scheme	Standard
	TX/RX delta frequency	0 Hz
	TX Nucleus	None
	TX delta frequency	0 Hz
	Coil elements	HE1-4
	Acquisition duration	0 ms
BOLD	GLM Statistics	Off
	Dynamic t-maps	Off
	Ignore meas. at start	0
	Ignore after transition	0
	Model transition states	Off
	Temp. highpass filter	Off
	Threshold	4.00
	Paradigm size	3
	Motion correction	Off
	Spatial filter	Off
	Delay in TR	0 ms
	Distortion Corr.	Off
	Contrasts	1